



Follow the instructions for complete installation of application.

- (1) First install DiabetesManager.msi by double clicking it.
After successful installation of above file, a folder 'DiabetesManager' will be created in Program files folder in drive C.
- (2) Now install redistributable package CRRedist2008_x86.msi from location C:\Program Files\DiabetesManager\ CRRedist2008_x86.msi
- (3) Shortcut DiabetesManager will be created on Desktop and in Program manager, if not then create a shortcut of exe file DiabetesManager.exe from Location "C:\Program Files\DiabetesManger"

Now your application is ready for operation.